

Standardise breast density analysis to provide personal care.

As awareness of breast density as a risk factor for cancer grows, so does the need for accurate analysis. Through a proprietary software algorithm powered by machine learning, the Quantra™ software analyses images for distribution and texture of parenchymal tissue and categorises breasts in four breast composition categories. This is critical for women with dense breasts as higher breast density is known to increase a woman's risk for breast cancer.

Quantra™ analyses both 2D and tomosynthesis images to help diagnose with confidence your most challenging breast images.




Improve assessment with unbiased algorithms powered by machine learning that analyse each patient's breast texture and pattern.



Elevate the standard of care and standardise reporting across the whole radiology practice.



View density on the AWS* to facilitate patient management protocols for adjunctive screening while the patient is in the office.




Example of a density score of "b" on the 3Dimensions™ Hologic acquisition workstation.

Scores are based on ACR BI-RADS Categories

Quantra™ software uses texture and pattern analysis for more consistent, more reliable scoring. It delivers the accurate information you need to confidently design patient-specific screening pathways and provide the highest quality care. There is a growing body of evidence that pattern and texture of fibroglandular tissue plays an equally important role as volume in mammographic cancer risk prediction.¹⁻³ The guidance language provided by the American College of Radiation (ACR) BI-RADS Atlas 5th⁴ Edition was changed to account for pattern and texture vs volume when determining density.

Quantra™ software:

- Overcomes subjectivity in visual assessment, which may impact patient care decisions.
- Produces an unbiased breast density assessment to support compliance with state specific density notification laws.
- Objectively supports managements of breast density protocols for adjunctive screening.
- Assigns breast density category based on analysis of parenchymal tissue distribution and texture
- Displays breast density scores on the Hologic acquisition workstation at the end of an exam.*
- Technologists can facilitate the implementation of high-risk/density-based patient management protocol.



Ordering details

Part Number	Description
DIM-LIC-QT	Quantra (2.2) license on AWS
DIM-LIC-QT-XFR	Quantra (2.2) license transfer to AWS
DIM-LIC-QT-UP-SC	Upgrade one Quantra license to 2.2 and transfer to AWS, service contract
DIM-LIC-QT-UP-NOSC	Upgrade one Quantra license to 2.2 and transfer to AWS, no service contract
DIM-LIC-QT-TRIAL	Quantra (2.2) 60 days demo license on AWS

*AWS with minimum 3Dimensions 2.1 or Dimensions 1.10 software

References:

1. Wang C, Brentnall AR, Cuzick J, et al. A novel and fully automated mammographic texture analysis for risk prediction: results from two case-control studies. Breast Cancer Res. 2017;19(1):114.
2. Winkel, RR, et al. 2016, BMC Cancer, Vol. 16, p. 414. - Mammographic density and structural features can individually and jointly contribute to breast cancer risk assessment in mammography screening: a case-control study.
3. Nielsen M, Vachon CM, Scott CG, Chernoff K, Karemore G, Karssemeijer N, Lillholm M, Karsdal MA. 2014, Breast Cancer Research, Vol. 16(2), p. R37. - Mammographic texture resemblance generalizes as an independent risk factor for breast cancer.
4. ACR BI-RADS Atlas 5th Edition. 2013. Available at: <https://www.acr.org/Clinical-Resources/Reporting-and-Data-Systems/Bi-Rads> (Accessed on 30/12/2019).

SS-00528-EUR-EN Rev.002 (03/20) Hologic Inc. ©2020 All rights reserved. Hologic, 3Dimensions, Dimensions, Quantra, Selenia, The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or international markets. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

www.hologic.com | info@hologic.com | +32.2.711.4680

No Compromise. No Comparison.