

Affirm[®] Prone Filler Pad

Transform your biopsy suite into a multi-use space with one simple pad

Efficient utilisation

For maximum flexibility in breast exams and biopsy procedures, the Affirm Prone Filler Pad provides a versatile multi-modality surface. This solid center pad converts the tabletop of the Affirm Prone Breast Biopsy System into a surface that can be used for breast biopsies or ultrasound screening and diagnostic exams.

The addition of an Affirm Prone Filler Pad transforms your biopsy suite into a versatile space that allows you to use the Affirm Prone table for Ultrasound when it is not being used for stereotactic biopsies.


Patient comfort

With the addition of the Filler Pad, the Affirm Prone Biopsy System can also be used for post-patient care in the supine position following ultrasound or stereotactic procedures. The memory foam padding is designed to optimise patient comfort.

Specifications

Width: 66 cm
Length: 91 cm
Depth: 4.1 cm
Weight: 1.8 kg

Note: The solid centre pad is only intended for use in conjunction with the 22.86 cm diameter table insert.



Contact your Hologic representative to learn more.

Breast & Skeletal Solutions | hologic.com | euinfo@hologic.com

Affirm[®]
Prone Biopsy System