

# TruNode® single-use gamma probe

**With the TruNode system, the single-use, wireless gamma probe from Hologic, you know with certainty.**

The TruNode system is the first high performance, sterile, single-use, wireless gamma probe that has been optimised for the most common radio-guided surgical procedures. The TruNode system is effective, in applications such as sentinel lymph node mapping and radioactive seed localisation.



## TruNode system features:

### Industry-Leading Resolution

The TruNode system features a 10% energy resolution, enabling the user to detect faint nodes closer to the injection site and facilitating precise hotspot localisation.

### Fingertip Control

Toggle between three modes of detection (high-sensitivity Scan, high-resolution Point and I-125 Seed) with a simple to use, fingertip controls from the sterile field.

### Three Detection Modes

With its three detection modes, the TruNode system optimises sentinel node, parathyroid and tumor localisation procedures.

### Minimal Incision Size

With a 10mm operative diameter, the TruNode system is designed to be minimally invasive and easy to maneuver.

### Sterile

The TruNode probe is fully sterile, eliminating drapes and reprocessing in the OR.

### Voice Reporting

The TruNode system's compact User Feedback Unit provides both audible and visual user feedback, reducing the need to monitor a display.

### Simple Placement

The TruNode system has a minimal footprint as the User Feedback Unit clips to any IV pole in the OR.

### Single-Patient Use

The TruNode system was designed to reduce the risk of cross-contamination and surgical site infections.

# A comprehensive **family of solutions**

The TruNode® single-use gamma probe fits right into Hologic's full range of breast conserving solutions. Hologic has a complete set of products for breast conserving surgery to help you improve efficiency, optimise the patient experience and deliver the best outcomes possible. Talk to your representative about the systems that can best suit your and your patients' needs, ensuring your patients have the best possible care when they need it most.

**BioZorb® 3D bioabsorbable marker**



**LOCalizer® wire-free guidance system**



Be sure to talk with your Hologic representative about which breast conserving solution is right for you and your patients.



Hologic BV, Da Vincilaan 5, 1930 Zaventem, Belgium.

Breast & Skeletal Health | [Hologic.com](https://www.hologic.com) | [euinfo@hologic.com](mailto:euinfo@hologic.com)

## TruNode®

Single-use Gamma Probe

MISC-05510-EUR-EN Rev 003 ©2022 Hologic Inc. All rights reserved. Hologic, The Science of Sure, BioZorb, LOCalizer, TruNode and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to [euinfo@hologic.com](mailto:euinfo@hologic.com).